

# Real-Time Drug Response Analysis of A549 Cells to Cisplatin with the CIMS-32 System

## Purpose

The purpose of this experiment is to evaluate the dose-dependent cytotoxic effects of cisplatin on A549 cells by monitoring impedance changes and determining the  $IC_{50}$  value using the CIMS-32 system.


## Experiment Overview

<b>Cell Name</b>	A549 cell line	<b>Seeding quantity</b>	3,000 cells/well
<b>System</b>	CIMS-32	<b>Cell chip</b>	CITO-16W01E-SGL
<b>Coating</b>	Poly-L-Lysine (PLL)	<b>Treatment</b>	Cisplatin

## Workflow

DIV0	DIV1-DIV3
<ul style="list-style-type: none"> <li>- Coat the CFPS plate</li> <li>- After seeding the cells, incubate for 1 hour and then add the media</li> </ul>	<ul style="list-style-type: none"> <li>- Impedance recording</li> <li>- Treat drug and record drug effects</li> </ul>

## Results



## Cisplatin-Induced Cytotoxicity in A549 Cells Monitored by Impedance Measurements

(A) Normalized impedance magnitude was measured at 17 kHz over time to evaluate the dose-dependent cytotoxic effects of cisplatin in A549 cells. Treatment was initiated at 24 hours, as indicated by the arrow, with concentrations ranging from 0.001 to 100  $\mu$ M. (B) The



maximum normalized impedance magnitude decreased in a dose-dependent manner, reflecting the cytotoxic effects of cisplatin. (C) Cell viability, determined using the MTT assay, was plotted against cisplatin concentration, yielding an  $IC_{50}$  value of 10.17  $\mu M$ . The  $IC_{50}$  value, calculated from impedance measurements, was 7.73  $\mu M$ , demonstrating close agreement with the  $IC_{50}$  value and confirming the consistent and reliable cytotoxicity measurements between the two methods. (D) Microscopic imaging further validated these findings by revealing dose-dependent morphological changes in A549 cells, consistent with the cytotoxic effects observed in both assays.

## Summary

- ✓ Cisplatin treatment induced dose-dependent cytotoxicity in A549 cells, as demonstrated by reduced impedance measurements with a  $IC_{50}$  value of 7.73  $\mu M$ , an  $IC_{50}$  value of 10.17  $\mu M$  determined by the MTT assay, and corresponding morphological changes observed microscopically. Compared to the MTT assay, the impedance-based method offers the added advantage of non-destructive, real-time monitoring. This measurement technique enables continuous observation of time-dependent cytotoxic effects and dynamic cellular responses to cisplatin treatment, effectively complementing the endpoint data provided by the MTT assay.

